

Introduction to adjusting and maintaining assistive technology

Medical Aids Subsidy Scheme

28 July 2022



Queensland
Government

Learning objectives

Explore common repairs and maintenance requirements

Orientate you to resources from the Medical Aids Subsidy Scheme and other sources about maintaining and adjusting equipment



This session's focus

Basic adjustable manual wheelchairs

Mobile shower commodes

Wheeled walking aids

Bathing and toileting aids

Pressure cushions

Sling hoists



Session structure

Evidence-based practice on repairs and maintenance

Sourcing instructional resources

Resources for manual wheelchairs, sling hoists, pressure cushions, mobile shower commodes and bath transfer benches

Maintenance checklists

Questions and discussions

Feedback form




How long does equipment last?

Equipment	Median replacement rate (according to equipment suppliers) in years ¹
Wheeled walking aid	5
Manual wheelchair	5
Manual wheelchair tyres / spokes	2-3
Hoist	5-9
Hoist batteries	2
Hoist sling	2.5
Pressure cushions	2.5-3.5
Mobile shower commode	5

1. Marini, I., Vang, C., Antol, D.L., et al. (2019). Redux: Empirical validation of medical equipment replacement schedules in life care plans. *Journal of Life Care Planning*, 17(2), 5-17.



Typical product performance



Wheelchair repair facts

44.8 - 88.3% full-time WC users completed a repair within 6 months^{1,4,5,6}

8.7% - 42% full-time WC users had an adverse consequence occur from wheelchair failure^{1,4,5,6}

Customisable MWC more likely to require repairs²

PWC more likely to require repairs^{1,2,3,4,5,7}

TIS more likely to require repairs^{2,3}

PWC users more likely to experience adverse consequences¹

Power seat functions not related to number of repairs¹

Mixed evidence whether power seat functions related to adverse consequences^{1,7}

Relationship between age / use of device and number of repairs is linear^{2,4}

Over time in the USA, number of repairs and adverse consequences is increasing⁷

1. McClure, L.A., Boninger, M.L., Oyster, M.L., et al. (2009). Wheelchair repairs, breakdown, and adverse consequences for people with traumatic spinal cord injury. *Archives of Physical Medicine & Rehabilitation*, 90(12), 2034-4.
2. James, A.M., Pramana, G., Schein, R.M., et al. (2022). A descriptive analysis of wheelchair repair registry data. *Assistive Technology*, online early.
3. Mhatre, A., Pearlman, J., Schmeler, M., Krider, B., & Fried, J. (2022). Community-based wheelchair caster failures call for improvements in quality and increased frequency of preventative maintenance. *Spinal Cord*, 60(1), 58-62.
4. Henderson, G.V., Boninger, M.L., Dicianno, B.E., & Worobey, L.A. (2022). Type and frequency of wheelchair repairs and resulting adverse consequences among veteran wheelchair users. *Disability & Rehabilitation: Assistive Technology*, 17(3), 331-337.
5. Worobey, L.A., Heinemann, A.W., Anderson, K.D., et al. (2022). Factors influencing incidence of wheelchair repairs and consequences among individuals with spinal cord injury. *Archives of Physical Medicine & Rehabilitation*, 103(4), 779-789.
6. Toro, M.L., Worobey, L., Boninger, M.L., Cooper, R.A., & Pearlman, J. (2016). Type and frequency of reported wheelchair repairs and related adverse consequences among people with spinal cord injury. *Archives of Physical Medicine & Rehabilitation*, 97(10), 1753-60.
7. Worobey, L., Oyster, M., Nemunatis, G., Cooper, R., & Boninger, M.L. (2012). Increases in wheelchair breakdowns, repairs, and adverse consequences for people with traumatic spinal cord injury. *American Journal of Physical Medicine & Rehabilitation*, 91(6), 463-469.

Common repairs by device^{1,2}

Category of device	Most common repairs (more to less frequent per category)
MWC	Wheels / castors / tires / forks General servicing Armrest Legrest / footrests
PWC	Batteries/cables Electronics General servicing Frames
Scooter	Batteries/cables General servicing

1. James, A.M., Pramana, G., Schein, R.M., et al. (2022). A descriptive analysis of wheelchair repair registry data. *Assistive Technology*, online early.
2. Toro, M.L., Worobey, L., Boninger, M.L., Cooper, R.A., & Pearlman, J. (2016). Type and frequency of reported wheelchair repairs and related adverse consequences among people with spinal cord injury. *Archives of Physical Medicine & Rehabilitation*, 97(10), 1753-60.



Consequences of wheelchair breakdowns^{1,2}

- Immobility
- Direct injury
- Worse pain, e.g., through increased rolling resistance with castor breakdown
- Poorer self-perceived health
- Higher risk of hospitalisation (odds ratio 2.17, $p < 0.10$)
- Higher risk of pressure injuries (odds ratio 1.72, $p < 0.05$)

1. Hogaboom, N.S., Worobey, L., Houlihan, B.V., Heinemann, A.W., & Boninger, M.L. (2018). Wheelchair breakdowns are associated with pain, pressure injuries, rehospitalization, and self-perceived health in full-time wheelchair users with spinal cord injury. *Archives of Physical Medicine & Rehabilitation*, 99(10), 1949-1956.

2. Wilson-Jene, J., Mhatre, A., Ott, J., et al. (2021). Rolling resistance of casters increases significantly after two years of simulated use. *Journal of Rehabilitation Assistive Technology Engineering*, 8, 20556683211025149.



How important is maintenance?

Abstract

Objective: To investigate whether active intervention using a compiled checklist for wheelchair check-ups increases user satisfaction and/or decreases accidents, near accidents and pressure sores.

Design: A randomized controlled trial comparing active intervention versus standard intervention for prescribed, manually propelled wheelchairs.

Setting: Patients within primary health care of Bora's and Bollebygd municipalities, a mixed urban and rural population.

Subjects: Users of manually propelled wheelchairs over 16 years of age.

Interventions: The accident rate, extent of pressure sores, number and extent of repairs, reconditioning, adjustments as well as user satisfaction were measured initially and at one year. In the standard intervention, the user and carer were encouraged to initiate contact when necessary. In the active intervention, an occupational therapist performed a scheduled, thorough check-up of the wheelchair, following a compiled checklist for safety, comfort and positioning, manoeuvrability and transportation.

Results: Of 253 registered wheelchair users, 216 were suitable and randomized. In the active intervention group, 99% (95% confidence interval 96 / 100%) of the inspected wheelchairs required maintenance. The incidence of accidents was unchanged in the standard intervention group, but decreased to 'no accidents' in the active intervention group ($p / 0.03$). User satisfaction was not affected by the active intervention.

Conclusion: Most wheelchair users are unable to determine on their own when adjustments are needed. An active check-up on manually propelled wheelchairs seems to reduce accidents. More information is available at <http://www.wheelchair.se>


- RCT manual wheelchair regular inspection (1 x 1hr per year) vs standard reporting
- Reduction in accidents with regular inspection
- Increased number of minor repairs by OT and exchange of parts



Toolkit

- Tape measure
- Socket wrench set
- Allen / hex key sets (metric and imperial)
- Adjustable spanners (small, medium, large) / wrench set
- Pliers (needle nose)
- Spoke key
- Bike pump ± tire pressure guage
- Quick unpick
- Screw drivers (Phillips and flat head)
- Lanolin-based lubricant
- Fabric protectant spray
- Upholstery protectant spray
- Puncture repair kit
- Velcro (adhesive)
- Scissors
- File
- Utility knife





MASS website resources: Mobility aid repairs and maintenance

Generic maintenance

Fact sheets

- [Tyre repair and replacement](#) (PDF, 182kB)

Video

- [Insurance](#) (WMV, 1.4MB)
- [Tyre maintenance](#) (WMV, 6.5MB)
- [Cleaning and lubrication](#) (WMV, 5.9MB)
- [Maintenance checks for electrical equipment](#) (WMV, 2.2MB)
- [Mechanical maintenance check - Part 1](#) (WMV, 5.5MB)
- [Mechanical maintenance check - Part 2](#) (WMV, 4.7MB)

Manual wheelchairs

Fact sheets

- [Adjusting your manual wheelchair, power wheelchair or mobile shower commode](#) (PDF, 761kB)

How to use a power wheelchair (PWC)

Fact sheets

- [You and Your Power Wheelchair - a guide to registration, insurance, maintenance and repair](#) (PDF, 502kB)
- [Adjusting your manual wheelchair, powerchair or mobile shower commode](#) (PDF, 761kB)

Video

- [Delivery and set up](#) (WMV, 2.1MB)
- [Battery charging](#) (WMV, 2.1MB)
- [Registration and insurance](#) (WMV, 3.2MB)
- [Use of the chair in wet or hot weather, indoor, outdoors](#) (WMV, 6.0MB)
- [Programming for individual needs](#) (WMV, 1.8MB)



MASS website resources: Daily living aid repairs maintenance

Hoist and sling use

Video

- [Hoist and sling use, maintenance checks \(WMV 6101 kB\)](#)
- [Battery charging and emergency lowering \(WMV 813 kB\)](#)

Mobile shower commodes

- [Adjusting your manual wheelchair, powerchair or mobile shower commode \(PDF 782 kB\)](#)

How to use pressure redistribution cushions, mattresses or overlays (PRC or PRM)

Video

- [Correct use of pressure redistribution equipment and skin checks \(WMV 8405 kB\)](#)
- [Offloading pressure \(WMV 4795 kB\)](#)
- [Alternating pressure mattresses and overlays \(WMV 1943 kB\)](#)
- [Cushion cleaning and maintenance \(WMV 2616 kB\)](#)
- [Mattress and mattress overlay cleaning, maintenance and faults \(WMV 2616 kB\)](#)




New MASS resources (SCU student project)



in the fork. If the wheel moves sideways or wobbles the bearings need service or replacement.

Manual Wheelchair Maintenance and Cleaning

[Manual Wheelchair Maintenance and Cleaning - YouTube](#)



Hoist Maintenance and Cleaning

[Hoist Maintenance and Cleaning - YouTube](#)

Completed by Mimi McCauley & Monique Nguyen




MASS workshop resources

Manual wheelchair - basic adjustments (Aspire Evoke 2)

Activity: Pick one of the clinical scenarios and adjust this manual wheelchair to be suitable for this client presentation.

Wheelchairs - what parts / many can be adjusted. user manual (online) for detailed adjustments to your chair.



Example clinical scenarios;

Vivian uses her wheelchair around a flat home. She has problems with balance in standing. She propels the wheelchair mainly with her feet. She often propels her wheelchair backrest and likes a firm backrest so she can push against it. She has arthritis in her wrists and has trouble with very tight brakes. She requires a 16" seat depth. Hint: low front and back seat heights, flat seat, upright backrest, medium tension brake, medium height backrest, medium axle position.

- Resources for inspecting, maintaining, adjusting, and cleaning a range of resources.
- Short videos
- Case scenarios

Try one
Scan here →



Sourcing other information on adjusting and maintaining

Online user manual

9. ADJUSTMENTS

Aspire Evoke 2 & Evoke 2 JNR Manual Wheelchairs have multiple adjustment points to ensure maximum comfort for users. It is very important that the below steps are followed when adjusting the wheelchair to ensure maximum safety.

Adjusting Backrest Height

1. Peel back upholstery in the lower edge of the backrest cane to reveal two bolt heads.
2. Remove 8 bolts (4 x each side) and reposition canes as required.
3. Tighten bolts and reattach upholstery.



Adjusting Backrest Tension

1. Peel back upholstery to reveal velcro adjustment straps.
2. Reposition straps to desired position and reapply upholstery securing in place with velcro.




Adjusting Seat Depth and Rear Wheel Position

1. Remove all bolts as shown and move backrest receiver bracket to the desired seat depth.
2. To adjust rear wheel position and height, the axle mount and receiver can be adjusted as shown. Once at the desired position, ensure all hardware is tightened securely.



NOTE: Ensure the rear wheels are located in symmetrical positions on both sides of the chair to maintain stability.

Online videos



How to adjust the Aquatec Ocean Ergo Shower Chair

1,500 views Dec 19, 2019

Like Dislike Share Download Save ...

InvacareCorp
5.28K subscribers

SUBSCRIBE


Sourcing other information on cleaning

On the item



e.g., Relax Duogel cushion

In the product manual



Online user manual

12. CARE & MAINTENANCE

Maintenance Schedule

Weekly	Periodically	Yearly	Serviceable Components
<ul style="list-style-type: none"> Check brakes are functioning properly. 	<ul style="list-style-type: none"> Check for signs of wear on front castors. Check for signs of wear on rear tyres. 	Service to be completed by authorised Aspire service agent	<ul style="list-style-type: none"> Wheels & Tyres Forks Axles Armrest receivers Seat sling Footrest mounts Footplates Armrest Pads

Care for Your Chair

- The seat and backrest upholstery can be removed by loosening the screws and wiping with a wet sponge and a mild detergent. All other components can be cleaned with a damp cloth.
- Check tyre condition at least twice per month.
- Check that rear wheel spokes are not loose. Loose wheel spokes can reduce the strength of the rims and must be safely adjusted by a qualified technician.
- Always keep rear wheel quick release axles clean and free from dirt and hair build up.
- Should push rims become scratched these should be de-burred or replaced to avoid damage to skin.
- Brakes should be checked regularly.

Should you have any doubts about the correct functioning of the wheelchair or of any of its components, please contact AIDACARE on 1300 133 120 or your local Aspire distributor.

The tool kit provided will allow basic adjustments to the wheelchair, for further information on repairs and service contact AIDACARE on 1300 133 120.


e.g., Aspire Ex

Maintaining & adjusting wheeled walkers

With your wheeled walking aid, you can;


- Change the handle height
- Adjust the brake tension
- Take off and clean the castor wheels

Depending on the walking aid you are using, it may differ from the following video.



[Wheeled walking aid adjustment](#)

brake tension. In the following video you will see a general introduction to inspecting and cleaning a 4-wheeled walker.




[Wheeled walking aid inspection and cleaning](#)


Maintaining manual wheelchairs



[Mobility aid resources](#)
[| MASS website](#)



[Manual Wheelchair Maintenance and Cleaning - YouTube](#)



[Wheelchair Service Training Package](#)

Adjusting manual wheelchairs

Adjustments to your manual wheelchair, power wheelchair or mobile shower commode

Preface
After taking delivery of a manual wheelchair, power wheelchair or mobile shower commode, several adjustments are often required to set the equipment up to fit the user and to obtain maximum benefit from use.

The following document uses a **suggested sequence of adjustments** starting with the footrest, then the seat and backrest, finishing with rear wheel position and balance of the chair. The sequence applies when considering the needs of the user.

FIRST set a seat height and angle for posture, transfers and shoulder to push rim relationship;

THEN move rear wheels forward / rearward to suit the stability needs of the user (Tomlinson 2000).

Often many of the seat adjustments have been decided when trialling the wheelchair, so these may be set up at the time of delivery by the supplier. Footrest height and back height often need adjustment at the time of delivery. Rear wheel position may be changed as the wheelchair user develops greater confidence or their needs change. Some **mechanical experience** is necessary to make most changes. The supplier, prescribing therapist or another person with suitable mechanical experience and access to tools may be able to assist.

The **product instruction manual** may have more detailed information on specific adjustments for equipment and should be referred to where available. This guide provides the rationale and aims for each adjustment, with instructions and product examples to give a general overview of adjustments for many different types of equipment.

A number of the adjustments discussed in this document will affect the overall stability of the wheelchair. As such it is important to carefully trial all changes to ensure the setup remains both functional and safe. In so doing consideration should be given to the environment of use with particular consideration of slopes and obstacles that are likely to be encountered.

Abbreviations Used	Manual wheelchair MWC	Power wheelchair PWC	Mobile shower commode MSC
Contents			
Adjustment	Aims	Applies to	Page
1. Footrest Height and Angle	<ul style="list-style-type: none"> Foot and lower limb support contributes to a stable pelvis, providing a base of support for the trunk Thighs, feet and buttocks all bearing weight Clearance under the footrest for outdoor mobility 	MWC, PWC, MSC	3
2. Seat Height and Angle	<ul style="list-style-type: none"> Provide clearance under the footrest and a seat height suited to propulsion, transfers and activity (tables, vehicle head height, etc). Seat angle contributes to stability in sitting 	MWC, PWC	4-7
3. Backrest Height and Seat to Back Angle	<ul style="list-style-type: none"> Balance the trunk over the pelvis. Not falling forward at the upper body. Pelvis does not slide forward. Promote normal spinal curves and prevent deformities To maximise reach and activity from the chair Trunk support decided after the pelvis and lower limbs are stable 	MWC, PWC	8-10
4. Rear Wheel Position (includes PWC seat position)	<ul style="list-style-type: none"> Make the chair easy to push and turn, allowing for stability needs and balance skills of user For PWCs the user's weight should be distributed to prevent tipping on slopes and for the drive system to perform as effectively as possible. For rear wheel drive chairs most user weight should be over the rear wheels aiming not to load the front castors too heavily, with mid or centre wheel drive chairs most user weight should be over the drive wheels 	MWC, some MSCs, some PWCs	11-12
5. MSC Seat Position	<ul style="list-style-type: none"> Move the seat so the aperture aligns with the toilet 	Some MSCs	13

Adjustments to your wheelchair, power wheelchair or mobile shower commode (Last reviewed June 2011) Page 1 of 13

Manual wheelchair - basic adjustments

Wheelchairs will vary greatly in what parts can be adjusted, and in many cases how they are adjusted. Refer to the product user manual (usually available online) for details about adjustments to your particular chair.

This activity uses the **Aspire Evoke 2 self-propelled wheelchair**. Some of the adjustments you can do include;

- Change the leg rest length
- Change the castor height (which changes the front seat height)
- Change the rear axle height (which changes the back seat

Manual wheelchair - basic adjustments (Ottobock Start M2S)

Activity: Pick one of the clinical scenarios and adjust this manual wheelchair to be suitable for this client presentation.

Wheelchairs will vary greatly in what parts can be adjusted, and in many cases how they are adjusted. Refer to the product user manual (usually available online) for details about adjustments to your particular chair.

WHEELCHAIR SERVICE TRAINING PACKAGE
Reference Manual for Participants
BASIC LEVEL

World Health Organization
USAID

PROBLEM: Wheelchair is too wide.
Sometimes the smallest wheelchair available is still too wide for the wheelchair user. If the wheelchair seat is too wide, the wheelchair user will find it difficult to sit upright, and is likely to collapse to one side. Some simple foam inserts can help to provide the wheelchair user with support to sit upright.
Solution: One solution is described below.

Solution	Caution
<ul style="list-style-type: none"> 1. Add foam inserts on each side of the pelvis. 2. Measure the space between the wheelchair user (sitting) seated in the wheelchair with the back against the backrest; use the side of each armrest. 3. Cut and add foam inserts to fill the space between each side of the wheelchair user and the side of the wheelchair. 4. Check fit. 5. Upholster and attach the foam inserts. If the armrests are solid, the inserts may be attached to the armrests. They may also be attached to the top of the cushion. 	<ul style="list-style-type: none"> The inserts only need to come as far forward as the user's trunk. Ensure the cushion matches the seat width. If the wheelchair is very wide, the wheelchair user will also have trouble reaching the push rims. Check for the push rims can cause a shoulder injury. Consider carefully whether the wheelchair is safe for the wheelchair user before proceeding with selecting it.

73

Mobility aid resources
MASS website

Manual wheelchair basic adjustment – Aspire Evoke 2

Manual wheelchair basic adjustment – Ottobock Start M2S

Wheelchair Service Training Package

Maintaining hoists



[Daily living resources | MASS Website](#)



Hoist Maintenance and Cleaning

[Hoist Maintenance and Cleaning - YouTube](#)

Maintenance Checklist – Hoists

Item	Task and Instructions
General	<ul style="list-style-type: none"><input type="checkbox"/> Check all parts for external damage or wear<input type="checkbox"/> Check all hardware, attachment points and stress parts, such as slings, hooks, and any pivot points for signs of wear, cracking, fraying, deformation, or deterioration<input type="checkbox"/> Verify that the hand control is functional (lifting and leg movements)<input type="checkbox"/> Check the emergency stop button
Batteries	<ul style="list-style-type: none"><input type="checkbox"/> Charge the battery regularly – if the hoist is used today, it should be charged today
Controls	<p>To ensure reliable operation of the hoist, check that:</p> <ul style="list-style-type: none"><input type="checkbox"/> All plugs are plugged in fully<input type="checkbox"/> Control or cord is not cracked, stretched, or damaged
Actuator and Spreader Bar	<p>Make a visual inspection to ensure:</p> <ul style="list-style-type: none"><input type="checkbox"/> Spreader bar is operating freely<input type="checkbox"/> High tensile pin attaching the spreader bar to the boom is moving freely but not excessively; not worn or rounded off at end<input type="checkbox"/> Bolt that secures the pin through the spreader bar to the boom is secure<input type="checkbox"/> Any clips, rings, or other means of attaching the actuator or spreader bar to the boom are intact and in good condition
Base legs, Mast, and Boom	<p>Make a visual inspection to ensure:</p> <ul style="list-style-type: none"><input type="checkbox"/> Legs can spread their full distance and lock in the open and closed positions<input type="checkbox"/> All bolts on the base legs are secure and none are missing<input type="checkbox"/> No visible cracks or bends in the base legs, mast, or boom


[Hoist & sling maintenance checklist | MASS](#)

Adjusting hoists

adjustment and troubleshooting

Activity: Try the following...

- Swap a standard yoke with a pivot frame yoke
- Try doing a transfer with each style of yoke
- Test the emergency stop button



Sling hoist adjustment (swapping a standard and pivot frame yoke) and troubleshooting

The image shows a person in a red shirt using a screwdriver to adjust a white hoist yoke. The yoke is part of a larger hoist assembly. The person is focused on the adjustment, and the background is a plain grey wall. The video player interface includes a play button in the center of the image and a small navigation icon in the bottom right corner.

Sling hoist adjustment (swapping a standard yoke and pivot frame yoke) and troubleshooting



Maintaining pressure cushions





[Daily living resources | MASS Website](#)

Cleaning and maintenance Checklist – Cushions

Please note:

- Cleaning and disinfecting are separate processes. Cleaning must precede disinfection.
- Before use by a different individual, clean, disinfect and check the product for proper functioning.
- If the integrity of the cushions have been compromised, do not use and seek replacement through MASS.

Inspection Checklist – Cleaning the cushions

Cushion	Task and Instructions
<p>'Action Gel' Action Pilot Gel Cushion</p> 	<p>Cushion cover:</p> <ul style="list-style-type: none"><input type="checkbox"/> Machine wash (60°) with Napisan <p>Cushion:</p> <ul style="list-style-type: none"><input type="checkbox"/> Gently scrub and wipe down with damp cloth and bleach or disinfectant (Viraclean) <p>Maintenance:</p> <ul style="list-style-type: none"><input type="checkbox"/> Check for any rips or damage to cushion cover<input type="checkbox"/> Check gel cushion for any large rips or tears in the gel
<p>'EQUAGEL' General Equagel Cushion</p> 	<p>Equagel cushion cover:</p> <ul style="list-style-type: none"><input type="checkbox"/> Cold machine wash with detergent<input type="checkbox"/> Air or tumble dry with NO heat <p>Equagel component:</p> <ul style="list-style-type: none"><input type="checkbox"/> For basic cleaning, hand wash in warm soapy water<input type="checkbox"/> For sanitation purposes, wash in a diluted bleach solution bath using 1 part household liquid bleach per 9 parts water <p>Maintenance:</p> <ul style="list-style-type: none"><input type="checkbox"/> Check pressure cushion and cover for any tears or damage

[Cushion Cleaning Checklist | MASS Action gel, Equagel, Relax Easy, ROHO](#)

Maintaining & cleaning pressure cushions¹



[Cleaning a ROHO Cushion - YouTube](#)



[Cleaning & maintaining a hybrid foam-gel pressure cushion \(Axiom SP\)](#)



[Action Products - YouTube](#)



[General advice for cleaning pressure cushions and covers](#)

Adjusting pressure cushions



[ROHO Cushions Set Up & Adjustment - YouTube](#)



[Adjusting a highly configurable foam pressure redistribution cushion \(Domino Inserto\)](#)



[Adjusting a Varilite cushion - YouTube](#)




Maintaining & adjusting bathing / toileting aids

maintaining a bath transfer bench

Activity:

- Watch the video below
- Check the bath transfer bench for wear and tear
- Set the seat height at the ideal height for you
- Swap the backrest over to the other side (e.g., if the client is facing the other way in the bath)



[Adjusting & maintaining a bath transfer bench](#)



Adjusting bathing / toileting aids

Adjusting a mobile shower commode - tilt-in-space K-care MSC

Activity:

- Complete the adjustments to the mobile shower commode as per video below
- Check, how easy is it to reach for perianal hygiene (either by self or a carer) with different seat configurations

Every mobile shower commode is different. Some will have a different range of adjustments possible. Some may use different tools (or allow tool less adjustment). In this case we'll be looking at a K-care tilt-in-space

[Adjusting a K-care tilt-in-space mobile shower commode](#)



[How to adjust the Aquatec Ocean Ergo Shower Chair - YouTube](#)

Maintenance checklists (peer-reviewed)

- Engineer cane/walker inspection
- Engineer wheelchair inspection
- Wheelchair maintenance assessment tool manual wheelchairs (p213-28)
- Wheelchair maintenance assessment tool power wheelchair (p238-54)
- Wheelchair Components Questionnaire¹
- Wheelchair maintenance checklist (p90-93)²

1. Rispin, K., Dittmer, M., McLean, J., & Wee, J. (2017). Preliminary reliability and internal consistency of the Wheelchair Components Questionnaire for Condition. *Disability & Rehabilitation: Assistive Technology*, 12(8), 852-856.
2. Sarcione, J.A., Kopec, C.L., & Trimby, A.J. (2004). Evaluating a wheelchair maintenance resource. Thesis, Worcester Polytechnic Institute.



Responsibilities & best practice

CLIENT

- Cleaning
- Simple maintenance
- Early detection and reporting of issues
- Proper use and following manufacturers instructions

PRESCRIBER

- Client education on care and maintenance
- Periodic equipment inspection and report/repair (ideally before warranty ends)
- Consider backup equipment / plans

MASS

- Arranging repairs and maintenance for wear and tear
- Providing replacement equipment at end of economically serviceable life
- Scheduled servicing of hoists and PWCs

Refer to the [MASS General Guidelines section 12](#)



Contacts

Dr Kieran Broome, Occupational Therapist

Advanced Clinical Educator, Medical Aids Subsidy Scheme

Bridget Manning, Speech Pathologist

Advanced Clinical Educator, Medical Aids Subsidy Scheme

Email: MASS-Education@health.qld.gov.au

Ph: 3136 3526

Local MASS service centre

For repairs

Brisbane email: MASS-Repairs@health.qld.gov.au

Brisbane phone: 3136 3525

Townsville email: MASS-Equipment-TSV@health.qld.gov.au

Townsville phone: 4433 8000



Feedback

Please complete our feedback form to receive a certificate of attendance.



[Link to feedback form](#)

