

Hollister Extended Wear Leg Bag

Continance Aids Product Information Sheet

Hollister Extended Wear Leg Bag – Sterile

(This is not a weekly change leg bag)



Description of Product: Hollister Extended Wear Leg Bag is supplied with adjustable extension tube which attaches directly to a sheath, indwelling catheter or night drainage bag. These bags are not designed to be changed weekly. They are designed to be changed only when the catheter is changed, i.e., every 8 weeks.

NB: Leg straps are supplied with these bags.

Available in:

Product Name	Product Code	Outlet Tap	Tubing Length	Total Capacity	Current MASS Distributor
Hollister Extended Wear Leg Bag	9348	Lever Tap	Adjustable	540mL	Independence
Hollister Extended Wear Leg Bag	9349	Lever Tap	Adjustable	900mL	Independence



© State of Queensland (Queensland Health) 2024
creativecommons.org/licenses/by/3.0/au

For more information contact:

Medical Aids Subsidy Scheme

PO Box 281

Cannon Hill QLD 4170 Australia

tel (07) 3136 3665

MASS-ContinenceAids@health.qld.gov.au

health.qld.gov.au/mass