

# The new Equipment Standing Offer Arrangement: new mobility and daily living aids available through MASS

Medical Aids Subsidy Scheme

31 July 2024



**Queensland**  
Government


# Learning objectives

- Understand the role of standing offer arrangements (SOA) in prescribing equipment
- Know features of new mobility and daily living aids on the new SOA
- Be aware of the new subsidy levels for mobility and daily living aids

# Transition period



# Development of new MASS Standing Offer Arrangement (SOA)



# Assistive Products not Included on the new SOA

Ceiling hoists

Specialised strollers

Sleep positioning systems

Bath hoists

Rigid frame manual wheelchairs

Multilift hoists

# New Temporary Guidelines

- New combined guidelines – all mobility aids and daily living aids are together: [Daily Living and Mobility Aids Temporary Guideline](#)
- Common statements can be found in the top section of each category, rather than spread throughout the document as previously e.g.

**For all bathing and toileting aids:**

Ensure Applicant/Carer is able to:

- Safely and effectively use the prescribed Assistive Product in all required areas of the home.
- Appropriately clean/maintain and store the prescribed Assistive Product.

- Notes at bottom of section that pertain to specific assistive products e.g.

- PWCs are provided for independent client use. Subsidy funding will not be approved for attendant propulsion.
- Basic PWC may be foldable/dismantlable for transport.

# Examples of Assistive Products no Longer on the Approved Product List

Ottobock  
wheelchairs

Oxford hoist

Reduced  
options for  
hoists

Arjo Sara  
Stedy

Kerry  
Multilift

Softform  
Premier  
Active 2

Juvo MSC  
(non height  
adjustable)

# Highlights of New SOA

- Large range of no quote sit and go manual wheelchairs
- Large range of basic power wheelchairs
- Seat elevate function on some scripted power wheelchairs
- Cost-effective options in most categories

## New, temporary subsidy funding

One year trial of tiered subsidy funding and then review

New subsidy levels reflect the current consumer group at MASS


Accessories are grouped together e.g. lateral support and headrest =  
postural supports

Clinical justification needed when requesting additional funding

Prescribers to contact MASS if their client falls outside of subsidy funding e.g.  
super-bariatric

With a co-payment of  $\leq$  \$200, an administrative waiver applies

# Tiered Approach to Subsidy Funding




## GOLDEN RULES

- Consider items on the SOA first (within the relevant category)
- Choose the most basic assistive product that meets your client's needs



# Trials of Assistive Equipment

- All equipment should be adequately trialled in the home environment, or it must be demonstrated how the equipment will be compatible with the home
- Trial is not required for the following products, where suitability can be confirmed with assessment of client function and review of specifications, with simulated trial or trial of a similar product where appropriate:
  - Bath Transfer Benches
  - Non-Standard Bathboards
  - Static shower chairs/Stools
  - Static and 3-in-1 Commodes and Over Toilet Aids
  - Wheeled Walking Aids
  - Basic Manual Wheelchairs – no accessories except anti-tips



- Post-delivery follow-up is to be provided by the prescriber or another AHP

# Remote Trial Subsidy

- Available for specified district postcodes
- Previous payment to supplier - \$200, now up to \$1000 (relates to distance from supplier)
- A range of products need to be taken for trial with the client
- A min. of two products per client needs to be trialled
- Every reasonable effort to explore alternative trialling options should be considered by the prescriber

# Quotes

- Moving away from needing a quote for each item e.g. no quote required for off-the-shelf products
- If not including a quote, need to ensure that product selection in MASS-eApply is accurate!
- A quote is required for all non-SOA and other scripted products
- Some manual wheelchairs may not require a quote



# Community Access

- Focus of MASS is to assist people to stay in their own homes
- However, there is a bit of flexibility with community access with wheeled walking aids

Lightweight frame may be considered to (one or more):


- Assist Applicant/Carer negotiating thresholds within the home.
- Be carried up/down stairs for use on different levels within the home.
- Assist with transportation for Applicant's medical appointments considering Applicant/Carer's existing health conditions/frailty.



# New Mobility and Daily Living Aids Product Spreadsheet

- One spreadsheet, not multiple, as in the past [MASS Equipment Approved Product List](#)
- Every model is listed, rather than one with different sizes
- Can use filters
- Two ways to locate cushions
  1. mobility aids (for a wheelchair)
  2. daily living aids (for a static chair)

➤ Screenshare





# Mobility Aids

Category – Mobility Aids	Primary Subsidy: Available towards products within the relevant category*	Primary Subsidy inclusions/Additional Funding: provide clinical justification to request consideration of funding*	Clinical Eligibility: Applicant has a permanent and stabilised condition or disability and requires Assistive Products to improve function and safety in the home and:
Wheeled Walking Aid (WWA)	\$450	Primary Subsidy includes frame, wheels, brakes.	Is unable to: Walk independently or effectively without support.
		Heavy Duty/non-standard sizing - \$250  Special Clinical Need - \$250	Requires a higher safe working load, stronger frame or non-standard sizing to <u>accommodate anthropometrics and/or usage.</u>  Has additional clinical need e.g. forearm/gutter support, lightweight frame. <u>Lightweight frame may be considered to (one or more):</u> <ul style="list-style-type: none"> <li>• Assist Applicant/Carer negotiating thresholds within the home.</li> <li>• Be carried up/down stairs for use on different levels within the home.</li> <li>• Assist with transportation for Applicant's medical appointments considering Applicant/Carer's existing health conditions/frailty.</li> </ul>

# Standard 4WW



## Staeligero

Aluminium  
7.5kg weight  
SWL 120kg  
51/55cm seat heights

## K-Care Economy 8"

Aluminium  
7.5kg weight  
SWL 136kg  
53cm seat height



## Hero walker

Aluminium  
6.7kg weight  
SWL 140kg  
55cm seat height

# Outdoor 4WW



## Freedom Stride

9kg weight  
SWL 135kg  
52cm seat height  
Side folding  
**11" front castor**



## Hero Outdoor Lite – Carbon Fibre

Standard/Tall  
7.1-7.3kg weight  
SWL 136kg  
52.5/60.5cm seat height  
Side folding  
**9.5" front castors**



## Hero Lite Portable Outdoors

8kg weight  
SWL 140kg  
55cm seat height  
Side folding  
**Flip up rear-wheels**  
**Integrated walking stick/crutch holder**

# Short/tall 4WW



Staebel Petit  
5.5kg weight  
SWL 175kg  
47cm seat height



Royale Premium 8"  
Low/Medium/Large  
8.6kg weight  
SWL 200kg  
50/55/59cm seat height



Hero Euro Lite Aluminium  
Low/Standard/Tall  
7.8-8.2kg weight  
SWL 130kg  
50/57/63cm seat height  
Side folding

# Adjustable seat height 4WW



Hero Ultra Light  
Adjustable height  
4.75kg weight  
SWL 113kg  
48-53cm seat height



Aspire Flex  
6/8" castors  
Adjustable height  
6.3-7.3kg weight  
SWL 136kg  
43-57cm seat height



Hero Deluxe with Adjustable Height  
Adjustable height  
6.7kg weight  
SWL 140kg  
46-59cm seat height

# Light-weight 4WW



K-Care Deluxe 6"  
Aluminium  
6.2kg weight  
SWL 160kg  
55cm seat height



Aspire Vogue Lightweight 2  
Aluminium  
6.45kg weight  
SWL 150kg  
54cm seat height



Rothcare Lightweight  
Carbon Fibre  
5.9kg weight  
Side folding



Hero Euro Lite– Carbon Fibre  
Low/Standard/Tall  
Side folding  
6.6-6.8kg weight  
SWL 130kg  
50/57/63cm seat height



Aspire Vogue Carbon Fibre  
Short/Medium/Tall  
Side folding  
6.2-6.35kg weight  
SWL 150kg  
50/55/62cm seat height

# Bariatric 4WW



K-Care Maxi 8"  
Seat width 580mm  
SWL 225kg



Royale Premium 8" -  
bariatric  
Seat width 550mm  
SWL 230kg

# Bariatric / Outdoor 4WW



Aspire Vogue Adventure  
10kg weight  
SWL 200kg  
**54cm seat width**  
53cm seat height  
Side folding  
**11" front castor**



Hero HD Outdoor Lite –  
Carbon Fibre  
Seat width 580mm  
SWL 204kg  
Side folding  
**9.5" front castor**



Hero HD Outdoor Lite -  
Aluminium  
Seat width 580mm  
SWL 204kg  
Side folding  
**9.5" front castor**

# Forearm walkers



Staelbel Navigator



Freedom Forearm Walker HD  
SWL 220kg

# Narrow 4WW



Trust Care Let's Go Indoor Walker  
Narrow – 52cm width  
Foldable

Lightweight (may attract additional subsidy)	Heavy duty (may attract additional subsidy)
Hero Medical Seat Walker Euro Lite	Peak Care Might Mack Walker
Aspire Vogue Carbon Fibre	Staebel Server HD
Aspire Vogue Lightweight	Aspire Vogue Adventure
K-Care Deluxe 6"	Hero Medical Seat Walker Outdoor Lite
	K-Care Maxi 8"
	Royale Premium 8"

Where to start with the most basic wheeled walker

Manual Wheelchair (MWC) – Basic	\$1,000	Primary Subsidy includes frame, seat, backrest, armrests, footrests, tyres/wheels, brakes.	<p>Is unable to: effectively mobilise/walk independently or with assistance within the home environment.</p> <p>Meets the following (one or more):</p> <ul style="list-style-type: none"> <li>Does not require complex seating or postural supports.</li> <li>Wheelchair is expected to be <u>primarily used to transit between areas in the home</u> rather than for prolonged sitting.</li> </ul>
		No script/quote required – product comes in prescribed seat width/depth without any accessories other than anti-tips.	
		Anti-tips - \$200	
		Heavy Duty/non-standard sizing - \$1,000	Requires a higher safe working load, stronger frame or non-standard sizing to accommodate anthropometrics and/or usage.

<p><b>Manual Wheelchair (MWC) – Scripted</b></p>	<p>\$2,500</p>	<p>Primary Subsidy includes frame, seat or seat pan, backrest or back canes, armrests, footrests, tyres/wheels, brakes.</p> <p>Requires specific MWC features, seating or postural supports unavailable in basic models.</p>	<p>Is unable to: effectively mobilise/walk independently or with assistance within the home environment.</p> <p>Meets the following:</p> <ul style="list-style-type: none"> <li>Spends prolonged periods seated in the wheelchair in the home; and</li> <li>Requires scripted wheelchair set up accessories/modifications, postural supports, and/or complex seating for efficiency of propulsion, positioning and/or function.</li> </ul>
		<p>Stump Support - \$650 Anti-tips - \$200</p> <p>Heavy Duty/non-standard sizing - \$1,000</p> <p>Other prescribed accessories/modifications/postural supports - \$2,000</p> <p>Backrest, including bracketry - \$1,500</p> <p>Customised/specialised frame - \$1,000</p>	<p>Specify if needed, no further justification required.</p> <p>Requires a higher safe working load, stronger frame or non-standard sizing to accommodate anthropometrics and/or usage.</p> <p>Requires (one or more):</p> <ul style="list-style-type: none"> <li>Modifications for posture, positioning, pressure care, function and/or fit within the home environment.</li> <li>Accessories such as height adjustable armrests, footrest options, padding, and postural supports e.g. lateral support, headrest.</li> </ul> <p>Requires (one or more):</p> <ul style="list-style-type: none"> <li>Support for posture and position.</li> <li>Tension adjustment and/or additional postural support to sustain sitting.</li> </ul> <p>Meets the following:</p> <ul style="list-style-type: none"> <li>Requires significant customisation or setup for postural supports and/or independent propulsion; and</li> <li>Utilises MWC &gt; 80% of the day e.g. has complex neurological condition; is active wheelchair user – spinal cord injury, bilateral amputee/or non-prosthetic user.</li> </ul>
		<p>Specialised Drive Mechanism (e.g. one arm drive) - \$1,000</p>	<p>Requires: Specialised drive mechanism to self-propel MWC e.g. dual rim or lever one-arm-drive, capstan push-rims.</p> <p>Note this does not include power assist drive mechanisms which are not covered by MASS subsidy funding.</p>

# Updated/upgraded model MWC

Aspire Lite  
Transit



Aspire  
Transit 2  
Armrests  
(swing away,  
removeable)  
16.4-16.6kg  
SWL 140kg




Karma  
8000  
Transit



Karma  
Transit 2  
Half-folding  
back  
Armrests  
(adjustabl)  
13.5kg  
SWL 130kg



Breezy  
Basix 2




Breezy  
Basix 2  
with half  
folding  
back  
SP and  
AP  
versions



# Updated/upgraded model MWC

Aspire Vida  
SP / AP



Aspire Vida  
X SP /  
Quickie  
Vida X SP



Ki Catalyst 4



Ki Catalyst  
4QS



Ki Catalyst  
4E



Company: \_\_\_\_\_

Ship To: \_\_\_\_\_

Account Number: \_\_\_\_\_ Date: \_\_\_\_\_

Address: \_\_\_\_\_

Quote  PO Number: \_\_\_\_\_

City: \_\_\_\_\_ State: \_\_\_\_\_

RTS: \_\_\_\_\_

Zip: \_\_\_\_\_ Phone: \_\_\_\_\_

Marked For: \_\_\_\_\_

Email/Fax: \_\_\_\_\_

Catalyst 4 K0004 \_\_\_\_\_ Charge  
113 Kg weight capacity

**Standard Inclusions:**  
 Indescent Navy Frame Color  
 Standard Back Upholstery  
 Flip back armrests  
 6 x 1 Inch front casters  
 24" Rear wheels  
 Composite Footplates

### Frame

#### Seat Width / Seat Depth

Select box corresponding to desired Seat Width and Depth

		Seat Width			
		16"	18"	20"	
Seat Depth	16"	N/C	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	18"	N/C	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

### Backrest

#### Back Post / Height Range / Back Height

Select Back Height in row corresponding to desired Back Post.  
 Angle Adjustable Back Post not available with Flip Back or Swing Away Armrests.

Back Post	Price	Back Height			
Standard	N/C	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17.5	<input type="checkbox"/> 18.5	<input type="checkbox"/> 19.5

#### Footrest / Extension / Seat to Footrest Length

Select Seat to Footrest Length in row corresponding to desired Footrest and Extension. Footrest Length of at least 2.5" shorter than Front Seat Height recommended. Extension Tube Hanger measurements based on 70° Hanger Bend.

Footrest	Extension	Seat to Footrest Length									
Extension Tube Hangers (70°)	Short	<input type="checkbox"/> 11	<input type="checkbox"/> 11.5	<input type="checkbox"/> 12	<input type="checkbox"/> 12.5	<input type="checkbox"/> 13	<input type="checkbox"/> 13.5	<input type="checkbox"/> 14	<input type="checkbox"/> 14.5	<input type="checkbox"/> 15	
	Medium	<input type="checkbox"/> 13	<input type="checkbox"/> 13.5	<input type="checkbox"/> 14	<input type="checkbox"/> 14.5	<input type="checkbox"/> 15	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17	
	Long	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17	<input type="checkbox"/> 17.5	<input type="checkbox"/> 18	<input type="checkbox"/> 18.5	<input type="checkbox"/> 19	<input type="checkbox"/> 19.5	
Elevating Leg Rests	Short	<input type="checkbox"/> 12.5	<input type="checkbox"/> 13	<input type="checkbox"/> 13.5	<input type="checkbox"/> 14	<input type="checkbox"/> 14.5	<input type="checkbox"/> 15	<input type="checkbox"/> 16	<input type="checkbox"/> 16.5		
	Long	<input type="checkbox"/> 13.5	<input type="checkbox"/> 14	<input type="checkbox"/> 14.5	<input type="checkbox"/> 15	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17	<input type="checkbox"/> 17.5	<input type="checkbox"/> 18

### Residual Limb Support

Residual Limb Support E1020  
 Left \_\_\_\_\_ Charge  Right \_\_\_\_\_ Charge

### Front Seat Height and Casters

#### Caster Size / Front Frame Type / Front Seat Height

Select Front Seat Height in row corresponding to desired caster size.

Caster Size	Front Seat Height			
6"	<input type="checkbox"/> 18	<input type="checkbox"/> 18.5	<input type="checkbox"/> 19	<input type="checkbox"/> 19.5

### Rear Seat Height and Wheels

#### Wheel Size / Tire Profile / Rear Seat Height

Select Rear Seat Height in row corresponding to desired Wheel Size and Tire Profile. Some Seat Depths may have conflicts with certain wheels in certain configurations.

Wheel Size	Tire Profile	Rear Seat Height			
24"	Full	<input type="checkbox"/> 16	<input type="checkbox"/> 17	<input type="checkbox"/> 18	<input type="checkbox"/> 19

#### Handrim

##### Standard Handrims

- Aluminum Anodized \_\_\_\_\_ Std
- Omit Handrim \_\_\_\_\_ N/C

#### Axle Bracket

- Offset Axle Bracket \_\_\_\_\_ Std
- Reverse Mount Offset Axle Bracket \_\_\_\_\_ N/C  
Not available with Swing Away Armrests

#### Axle

- Fixed Axle \_\_\_\_\_ Std
- Quick Release Axle \_\_\_\_\_ N/C

#### Anti-Tip

- Rear Anti-Tips E0071 \_\_\_\_\_ Charge

#### Additional Information:

# Budget MWC



Days Swift  
Lightweight  
(12.9kg total,  
7.8kg frame)



Hero Super Budget  
SP/AP  
17.6/15.4kg



Staebel Thrifty  
Transit  
Weight ~16kg

# Standard basic MWC



Karma S-Ergo 125  
Curved seating



Hero Premium Comfort  
"Comfy" Lite  
SP/AP  
SP weight 11.7-12.1kg  
AP weight 9.7-10.1kg



Staebel Deluxe Alloy  
PU armrest pads  
Height adjustable armrests  
Flip-down back  
Optional: Headrest



Staebel Emperor  
Transit (AP)  
Optional: elevating  
legrests

# Lightweight MWC



Karma Ergo Lite  
Transit

Weight 8.7-8.8kg  
SWL 100kg  
**Larger rear castor**



Hero Featherlite  
AP

Ultra lightweight 8.4kg



Hero Active Ultra Lite  
AP

Carbon fibre  
**Larger rear castor**  
Lightweight 8.9kg



Staebel Liberation  
Transit (AP)

Ultralight 7kg  
**Highly collapsible**

# Bariatric / heavy duty MWC



Freedom HD Wide

Heavy Duty

Seat width 560mm

SWL 250kg



Quickie M6

Heavy Duty

Seat width 560-750mm

SWL 295kg



Staelbel Gusto

Heavy duty

Seat width 610-710mm

SWL 240-290kg

# Highly adjustable / configurable MWC



Meyra Eurochair  
Highly configurable  
(including depth  
adjustable armrests)



Staelbel Sovereign  
Highly adjustable  
Multiple widths  
Optional: headrest  
Optional: elevating legrest  
Optional: stump support

# Highly scripted MWC



Quickie QXi  
Highly configurable  
Scripted  
Quick adjust



Quickie 2 / Quickie 2 Lite  
Highly configurable  
Lightweight (from 11.5kg)

Lower cost basic	Lower cost scripted / configurable	Lower cost heavy duty
Hero Medical Super Budget (AP/SP)	Breezy Basix (AP/SP)	Freedom HD wide wheelchair (SP)
Hero Medical Featherlite (AP)	Evoke 2 (SP)	Evoke2 HD (SP)
Days Swift (AP/SP)	Staebel Sovereign (SP)	Staebel Deluxe Alloy Wheelchair (SP)
Aspire Transit 2 (AP)	Meyra Eurochair (SP)	Staebel Gusto (SP) – for extra wide / high SWL
Aspire Vida SP		

Where to start with the most basic MWC

<p><b>Tilt in Space Manual Wheelchair (TISMWC) – Basic</b></p>	<p>\$4,000</p>	<p>Primary Subsidy includes tilt in space frame, seat, backrest, armrests, footrests, tyres/wheels, brakes, tilt mechanism.</p> <p>No script/quote required – product comes in prescribed seat width/depth without any accessories other than anti-tips.</p>	<p>Is unable to: effectively mobilise/walk independently or with assistance within the home environment.</p> <p>Requires tilt in space to (one or more):</p> <ul style="list-style-type: none"> <li>• Achieve or maintain a more upright, stable supported sitting position when head and trunk control are impaired.</li> <li>• Redistribute pressure/reduce concentrated areas of soft tissue on the seat.</li> <li>• Facilitate optimal positioning for comfort and function due to deformity/pain/involuntary movement/abnormal tone/seizure activity.</li> <li>• Better manage respiration, gastrointestinal function and/or facilitate hoist transfers.</li> </ul>
<p><b>Tilt in Space Manual Wheelchair (TISMWC) – Scripted/Tilt in Space Specialised Stroller</b></p>	<p>\$6,000</p>	<p>\$6,000 Includes frame, seat or seat pan, backrest or back canes, armrests, footrests, tyres/wheels, brakes, tilt mechanism.</p> <p>Requires specific MWC features, seating or postural supports unavailable in basic models.</p>	<p>Is unable to: effectively mobilise/walk independently or with assistance within the home environment.</p> <p>Requires</p> <ul style="list-style-type: none"> <li>• Scripted wheelchair set up accessories/modifications, postural supports, and/or <u>complex seating for positioning and/or function</u>; and</li> <li>• Tilt in space to (one or more): <ul style="list-style-type: none"> <li>– Achieve or maintain a more upright, stable supported sitting position when head and trunk control are impaired.</li> <li>– Redistribute pressure/reduce concentrated areas of soft tissue on the seat.</li> <li>– Facilitate optimal positioning for comfort and function due to deformity/pain/involuntary movement/abnormal tone/seizure activity.</li> <li>– Better manage respiration, gastrointestinal function and/or facilitate hoist transfers.</li> </ul> </li> </ul>

- Plus additional subsidies for:
- Heavy duty
  - Stump support
  - Anti-tips
  - Backrests
  - Other accessories for posture, position & pressure care

# TIS MWC



Staelbel Relax Tilt & Recline

Gas struts  
Manual leg elevate



Netti 4U CE Plus

Gas struts  
Variable backs  
Adjustable seat width



Breezy Relax2

Highly adjustable  
(seat height & depth,  
back height,  
armpads)

# TIS MWC



Karma VIP2 Tilt & Recline Self Propel  
Foldable  
Lightweight 17.4kg



Karma VIP515 Tilt Self Propel  
Foldable  
Lightweight 19.3kg



Aspire Rehab Rx Tilt & Recline  
Adjustable / scripted  
(MASS AP only)

# TIS MWC



Quickie Iris  
Scripted  
Rotation-in-space



Quickie SR45  
Scripted  
Rotation-in-space



Ki Liberty QS  
Quick supply

Ki Liberty FT  
Highly adjustable  
Gus struts  
Folding  
Self-propelled  
Heaviest part 11.7kg

Lower cost basic	Lower cost configurable	Lower cost heavy duty (SWL > 150kg)
Staebel Relax Tilt & Recline	Aspire Rehab Rx Tilt & Recline	Aspire Rehab Rx Tilt & Recline
Netti 4U CE Plus	Karma VIP515	
	Breezy Relax2 TIS	
	Karma VIP2 Tilt & Recline	

Where to start with the most basic TIS MWC

Power Wheelchair (PWC) – Basic	\$4,500	Primary Subsidy includes seat, backrest, armrests, footrests, motors, brakes, wheels/tyres, battery, charger, joystick.	<p>Unable to: Effectively mobilise/walk independently or with assistance or self-propel a manual wheelchair within the home environment.</p> <p>Meets all of the following:</p> <ul style="list-style-type: none"> <li>• Demonstrated ability to develop the skills necessary to effectively control and manoeuvre the requested PWC in all required areas of the home.</li> <li>• Does not require complex seating, postural supports or features e.g., tilt in space, seat elevate.</li> <li>• Does not need to negotiate variable surfaces.</li> <li>• PWC is expected to be primarily used to transit between areas in the home rather than for prolonged sitting.</li> </ul>
		Heavy Duty/non-standard sizing - \$1,000	Requires a higher safe working load, stronger frame or non-standard sizing to accommodate anthropometrics and/or usage.
	<p>Note:</p> <ul style="list-style-type: none"> <li>• PWCs are provided for independent client use. Subsidy funding will not be approved for attendant propulsion.</li> <li>• Basic PWC may be foldable/dismantlable for transport.</li> </ul>		

**Power  
Wheelchair  
(PWC) – Scripted**

\$10,000

Primary Subsidy includes seat or seat pan, backrest or back canes, armrests, footrests, motors, wheels/tyres, brakes, battery, charger, joystick.

Is unable to: effectively mobilise/walk independently or with assistance or self-propel a manual wheelchair within the home environment.

Meets all the following:

- Demonstrated ability to develop the skills necessary to effectively control and manoeuvre the requested PWC in all required areas of the home.
- Spends prolonged periods seated in the wheelchair in the home.
- Requires scripted wheelchair set up, features e.g. tilt in space/seat elevate, accessories/modifications, postural supports, and/or complex seating for positioning and/or function.

May be required to: Negotiate variable surfaces.

Plus additional subsidies for:

- Heavy duty
- TIS
- Elevate
- Accessories for posture, positioning, pressure care, function & fit
- Backrest
- Specialised controls
- Specialised frames

# Standard basic PWC



Merits Maverick 10  
MWD

# Bariatric / heavy duty PWC



Pride Jazzy 1450

Heavy Duty  
Seat width up to  
81.2cm  
SWL 272.1kg



Quantum 1450

FWD  
Seat width up to  
81.2cm  
SWL 272kg



Permobil M300 HD

SWL 204kg  
Seat width: up to  
63.5cm

# PWC (with all terrain potential)



Magic 360

Compact

Wheels: Urban /  
Crossover / Off-road



Magic Frontier V4

FWD/RWD

Wheels: Urban /  
Crossover / Off-road



Magic Frontier V6

MWD

Wheels: Urban /  
Crossover / Off-road

# PWC (with all terrain potential)



Quantum Edge Z

MWD

Off-road

(independent  
castors)



Quantum R-Trak

4pt independent  
suspension



Quantum 4Front 2

FWD

Suspension

Standard: 136kg SWL

HDL: 204kg SWL

Ant tilt avail

# PWC (with catapult function)



Dietz Sango  
Advanced  
FWD/MWD/RWD  
Option: catapult  
function



Meyra Orbits  
MWD  
Basic or complex  
Compact  
Options: Dual or single rear castor  
Optional: catapult

# Narrow PWC



Dietz Sango Slimline  
FWD/MWD/RWS  
Narrow (total width from  
53cm)  
Option: catapult function



Quantum Edge 3  
Stretto  
MWD  
Narrow (width from  
52.7cm)

# PWC (with unique feature/s)



Merits Dualer  
Hybrid FWD & RWD



TA iQ  
MWD/FWD/RWD  
Minimalist design  
Compact  
Low seat height (similar to MWC)  
Most power functions



Royale Travel Lite  
8/10/12" wheels  
Foldable  
Tilt/Recline  
~30kg with batteries removed

# Scripted PWC



Permobil M3  
MWD  
Ant tilt avail



Permobil F3  
FWD  
Ant tilt avail



Quantum Edge 3  
MWD



Quickie 500  
MWD/FWD/RWD  
Seating: Sedeo Pro/Lite  
Wheels: Urban /  
Crossover / Off-road

# Lightweight foldable PWC



Pride iGo

RWD

Foldable

19.8kgs with batteries  
removed



Pride iGo Carbon Fibre

RWD

Foldable

17.8kg with batteries  
removed



Pride Litestream –

Carbon Fibre

RWD

Foldable

16kg with batteries  
removed



Quickie Q50R

RWD

Carbon Fibre

Foldable

14.5kg with batteries  
removed

Basic PWCs within the primary subsidy	Lower cost scripted PWCs	Lower cost heavy duty PWCs (SWL > 150kg)
Pride Go Chair	Quickie Q300	Pride Jazzy Elite HD (SWL 204kg)
Pride iGo	Quickie Q400 with captain seat	Q700 with Sedeo Lite seating
Royale Medical Travel Lite	Quickie Q500 with sedeo lite seat	Magic Frontier V4/V6
Jazzy Select 6	Quickie Q700 with captain seat	Meyra MC2 Basic
Quickie Q50R Carbon Fibre	Permobil M1	Pride Jazzy 1450 (SWL 272kg)
Quickie Q100		
Pride Litestream		

Where to start with the most basic PWC



# **Pressure Redistribution Support Surfaces**

**Pressure Redistribution Cushion (PRC) – Wheelchair**

\$750

Primary subsidy includes standard cover, pump where relevant, all components required for operation of product.

Eligible for a manual or power wheelchair and meets the following (one or more):

- At risk of developing a pressure injury and/or is unable to effectively redistribute pressure.
- Current/history of pressure injury on areas in contact with the seat, determined not IAD – Incontinence Acquired Dermatitis.
- Major fixed skeletal deformity and/or motor/sensory loss with potential for pressure injury development.

Special Clinical Need - \$450

Requires customised or specialised cushion e.g. motor/sensory loss with potential for pressure injury development, higher safe working load or non-standard sizing to accommodate anthropometrics.

Prescribed accessories/modifications - \$500 e.g. wedges, fluid inserts, cutouts, incontinence cover.

Requires (one or more):

- Additional postural or positioning supports to accommodate/correct skeletal deformity and/or postural asymmetry.
- Specialised cover or inserts to meet continence/wound needs.

**Pressure Redistribution Cushion (PRC) – Static Chair**

\$750

Primary Subsidy includes standard cover, pump where relevant, all components required for operation of product.

Meets the following (one or more):

- At risk of developing a pressure injury and/or is unable to effectively redistribute pressure.
- Current/history of pressure injury on areas in contact with the seat, determined not IAD – Incontinence Acquired Dermatitis)
- Major fixed skeletal deformity and/or motor/sensory loss with potential for pressure injury development.

Special Clinical Need - \$450

Requires customised or specialised cushion e.g. motor/sensory loss with potential for pressure injury development, higher safe working load or non-standard sizing to accommodate anthropometrics.

Prescribed accessories/modifications - \$500 e.g. wedges, fluid inserts, cutouts, incontinence cover.

Requires (one or more):

- Additional postural or positioning supports to accommodate/correct skeletal deformity and/or postural asymmetry.
- Specialised cover or inserts to meet continence/wound needs.

# PRC




Action  
Shearsmart



Action Xact



Blake Medical  
(silver / gold)



Blake Medical  
Hybrid



Ki Axiom  
Adjustable



Comfort  
Curve



Comfort  
Embrace



Comfort M2



Trulife  
Relax Gel  
Cell

# PRC



Jay 2



Jay 2 Deep Contour




Jay Basic



Jay Fusion



Jay Ion




Jay Union



Jay Low-zone Eco



Jay Low-zone Switch



Jay Low-zone Gel Switch



Spex Standard




Spex High Contour




Spex SuperHigh Contour


# PRC



Stealth  
Pride Tru  
Comfort




Stealth Pride  
Spectrum Gel



Stealth Pride  
Spectrum  
(foam)



Stealth  
Pride  
Essense



Stealth  
Pride  
Glacial



Stealth Pride  
Premiere



Stealth  
Pride  
Simplicity



Stealth  
Pride  
Solution

# PRC



ROHO  
Contour



ROHO  
Enhancer



ROHO Mosaic



ROHO  
Hybrid  
Select



Staebel  
Bariatric



Vicair  
Adjuster 10

<b>Pressure Redistribution Mattress/Overlay (PRM), or Sleep Positioning System (SPS)</b>	\$4,000	Primary Subsidy includes standard cover, pump where relevant, all components required for operation of product.	Meets the following (one or more): <ul style="list-style-type: none"> <li>• At risk of developing a pressure injury and/or is unable to effectively redistribute pressure.</li> <li>• Current/history of pressure injury on areas in contact with the mattress, determined not IAD – Incontinence Acquired Dermatitis).</li> <li>• Major fixed skeletal deformity and/or motor/sensory loss with potential for pressure injury development.</li> </ul> Sleep Positioning System additional criteria (meets one or more): <ul style="list-style-type: none"> <li>• Requires support and positioning to facilitate improved respiration and/or swallowing.</li> <li>• Prevention of pressure injury through specific positioning.</li> <li>• Prevention of contractures and/or deformities.</li> </ul>
		<u>Special Clinical Need - \$2,000</u>	Meets the following (one or more): <ul style="list-style-type: none"> <li>• Higher safe working load or king single size to accommodate anthropometrics, bed mobility and/or existing bed size.</li> <li>• Significant history of pressure injury e.g. post-surgical management of pressure injury.</li> <li>• Severe restriction in mobility e.g. high-level quadriplegia, locked in syndrome.</li> <li>• More basic/cost effective options have been trialled/considered and demonstrated to be unsuitable.</li> </ul>
	<b>Note:</b> <ul style="list-style-type: none"> <li>• MASS subsidy is available for single size or king single size only. Larger mattresses/overlays may be requested – MASS subsidy may be applied for the equivalent single/king single option. The applicant will be required to contribute a copayment for the difference.</li> <li>• Not available for pain management only.</li> <li>• Sleep Positioning Systems are intended for applicants with complex positioning needs e.g. neurological condition such as acquired brain injury, cerebral palsy.</li> <li>• Funding cannot be applied towards the purchase of a bed or wedges/heel supports.</li> </ul>		

# PRM - Foam




AreaCare  
Multilayer  
foam




AreaCare Optimum  
Multilayer castellated  
foam with heel slope



AreaCare Premium  
Multilayer foam  
with heel slope



AreaCare Plus  
Multilayer  
foam



AreaCare Pro  
Multilayer foam



My Comfort  
Waterproof  
Foam



My Comfort Pressure  
Plus Sloped Heel  
Foam

# PRM – air & alt air



Better Living Overlay  
1:2 alt air



Better Living Mattress  
1:2 alt air



Procair Trio  
1:3 alt air



ROHO Prodigy  
Static air overlay




ROHO Sofflex 2  
Static air overlay


# PRM - hybrids




Forte Zephair Hybrid  
Hybrid alt air – foam



Lifecomfort Aspire Hybrid  
Hybrid alt air - foam



Carilex Hybrid  
Hybrid alt air - foam



Forte Alaise  
Hybrid air - foam



CairMax Duo  
Hybrid air-foam



# Daily Living Aids

<b>Patient Transfer Platform</b>	\$2,500	Primary Subsidy includes frame	Is unable to (one or more): <ul style="list-style-type: none"><li>• Effectively complete standing or non-standing transfer using grab rails and/or transfer devices such as a slideboard.</li><li>• Effectively complete standing or non-standing transfer with assistance.</li></ul>
----------------------------------	---------	--------------------------------	---

# PTP



Oxford Up Stand-Aid  
SWL 200kg  
Base 633x960mm



Etac Turner Pro  
SWL 200kg  
Height 800-1200mm



Molift Raiser Pro  
SWL 170kg  
Height 1200mm

Standing Hoist	\$4,100	Primary Subsidy includes frame, battery, charger, standard spreader, controller.	Is unable to (one or more): <ul style="list-style-type: none"> <li>Effectively complete standing or non-standing transfer using grab rails and/or transfer devices such as a slideboard.</li> <li>Effectively complete standing or non-standing transfer with assistance.</li> </ul>
Mobile Floor Hoist	\$3,500	Primary Subsidy includes frame, battery, charger, standard spreader, controller.  Non-standard attachment – pivot/power pivot or 4-point spreader - \$700	Is unable to (one or more): <ul style="list-style-type: none"> <li>Effectively complete standing or non-standing transfer using grab rails and/or transfer devices such as a slideboard.</li> <li>Effectively complete standing or non-standing transfer with assistance.</li> </ul> <u>Standard spreader is unsuitable.</u> Non-Standard attachment is required for comfort and/or positioning between lying and sitting and between surfaces e.g. bed to wheelchair.
		Other prescribed requirements e.g. increased lift, reach, safe working load - \$2,000	Requires a higher safe working load, increased lift range/reach to accommodate anthropometrics and/or complete transfers. E.g. reach and leg spread for transfer to recliner.
	Multilift Hoist	\$5,000	Primary Subsidy includes frame, battery, charger, standard spreader, controller.  Non-standard attachment – pivot/power pivot or 4-point spreader - \$700

Note: The function of the full lift component of the multi-lift hoist must be considered for current and likely future needs.

# Hoists



ISA Stand up lifter  
SWL 180kg



Birdie Evo Plus  
Electric leg  
SWL 180kg



Birdie Evo Xplus  
Electric leg  
SWL 210kg



Molift Mover 180  
SWL 180kg

<b>Mobile Shower Commode (MSC) Basic/Standard</b>	\$2,000	Primary Subsidy includes frame, <u>basic handmade/ seamless padded seat, backrest, armrests, footrests, wheels/tyres, pan and carrier, brakes and/or locking casters.</u>	Unable to: Effectively walk to and/or transfer to the toilet and/or shower.
<b>Tilt in Space Mobile Shower Commode (TIS MSC)/ Shower Trolley</b>	\$3,500	Primary Subsidy includes tilt in space frame, basic handmade/seamless padded seat, backrest, armrests, footrests, wheels/tyres, pan and carrier, brakes and/or locking casters.	Is unable to: Effectively walk to and/or transfer to the toilet and/or shower. Requires tilt in space to (one or more): <ul style="list-style-type: none"> <li>• Achieve or maintain a suitable posture for hygiene.</li> <li>• Redistribute pressure/reduce concentrated areas of soft tissue on the seat.</li> <li>• Better manage respiration, gastrointestinal function and/or facilitate hoist transfers.</li> </ul>

- Plus additional subsidies for:
- Height adjustable
  - Self-propulsion
  - Heavy duty
  - Stump support
  - Anti-tips
  - Backrests
  - Other accessories for posture, position & pressure care

<u>Height Adjustable \$500</u>	Non-standard height setting is required to facilitate transfers, positioning, or fit over toilet.
<u>Self-Propelled - \$500</u>	Large rear wheels required for self-propulsion or assisted access e.g. to traverse thresholds/small changes in level.

# MSC



Juvo Folding Shower  
Commode  
AP/SP  
Dismantleable



Juvo Tilt & Recline  
Shower Chair  
SWL 200/300kg



Better Living MSC  
Weight 12kg  
SWL 136kg



Aspire Aluminium  
AP/SP  
AP Weight 15.1kg  
SP Weight 20.5kg  
SWL 136kg  
Dismantleable

# MSC



Raz Z333 TIS  
Orbital tilt



Raz EZPZ100  
AP



Raz EZPZ200  
SP



Raz EZPZ300  
TIS  
Hinged tilt

Lower cost	Lower cost TIS	Lower cost heavy duty (SWL > 150kg)
Aspire aluminium (AP/SP)	Aquatec Ocean Ergo VIP	4HC Verto – height adjustable (AP)
Better Living (AP)	Raz EZPZ300	Aspire (AP/SP)
4HC Verto – height adjustable (AP)	Juvo Tilting	KCare Transit (AP)
Aquatec Ocean Ergo Transit (AP)	Aspire Tilt & Recline (AP)	Juvo (AP/SP)

Where to start with the most basic MSC

<b>Bath Transfer Bench, or Swivel Bath Seat, or Bath Lift, or</b>	\$500	Primary Subsidy Includes frame, seat, non-slip feet (if applicable).	Is unable to: Effectively transfer into a bath or shower and/or stand to shower.
		Suction feet per pair - \$70 Leg extensions per pairs - \$70	Specify if needed, no further justification required.

# BTB



Royale Medical  
Non-Padded  
Premium  
Aluminium  
SWL 280kg



Royale Medical  
Sliding & Swivel  
Aluminium  
SWL 130kg

# Swivel bathseat



## Hero Medical Swivel Bather

Overall width 710mm

Seat width 440mm

SWL 113kg

<b>Static Commode or Static 3 in 1 Commode or Static Over Toilet Aid</b>	\$250	Primary Subsidy includes frame, seat, non-slip feet.	Is unable to: Effectively walk to and/or transfer on/off the toilet in the home.  Note: The item may be required for bedroom or close-proximity toilet use or positioned over the toilet. A static 3 in 1 commode may also be positioned in the shower.
		Heavy Duty/non-standard sizing - \$250.	Requires a higher safe working load or non-standard sizing to accommodate anthropometrics.
		<u>Special Clinical Need - \$250</u>	Has additional clinical need e.g. padded seat, backrest.

# Commodes



Staebel Basic Bedside  
Commode with Padded  
Armrest

Seat width 440mm  
SWL 125kg



Hero Medical Bedside  
Commode

Seat width 440mm  
SWL 150kg



Aspire Maxi Bedside  
Commode

Seat width 600/650mm  
SWL 210kg

# 3-in-1 commode



## Juvo Trix

New sizes – seat width  
46/51cm  
SWL400kg



## Hero Medical Bariatric

(heavy duty all-in-one  
with padded back)  
Seat width 60cm  
SWL 225kg



## Hero Medical Drop Arm

### Bariatric

Seat width 65cm  
SWL 225kg

# OTA



Freedom OTA  
HBA370

Seat width 480mm  
SWL 160kg



Freedom OTA  
HBA400

Seat width 500mm  
SWL 140kg



Freedom OTA  
HBA406

Seat width 460mm  
SWL 200kg

# OTA



Juvo/Staebel Extracare OTF  
500/560/610/660mm seat width  
SWL 150/175/200/200kg



Juvo/Staebel OTF  
Aluminium



Etac Swift  
Seat width 540mm  
SWL 160kg

**Static Shower  
Chair/Stool**

\$250

Primary Subsidy includes frame,  
seat, non-slip feet.

Is unable to: Effectively and safely stand to shower.

Heavy Duty/non-standard sizing -  
\$250

Requires a higher safe working load or non-standard sizing to  
accommodate anthropometrics.

Special Clinical Need - \$250

Has additional clinical need e.g. padded/swivel seat, backrest.

# Shower stools / chairs



Etac Edge  
SWL 130kg



Etac Easy  
SWL 200kg



Etac Smart  
SWL 200kg



Better Living  
Shower Chair  
Aluminium  
Seat width 445mm  
SWL 110kg



Better Living  
Shower Stool  
Aluminium  
Seat width  
370mm  
SWL 120kg

# Shower stools / chairs



Freedom Shower  
Chair  
HBA402  
Seat width 490mm  
SWL 140kg



Freedom Shower  
Chair  
HBA407  
Seat width 460mm  
SWL 200kg



Freedom Shower Stool  
HBA403  
Seat width 490mm  
SWL 140kg





Freedom Shower Stool  
HBA408  
Seat width 400mm  
SWL 200kg

# New SOA - Application Process

MASS-eApply only

For more information and to register, visit:  
<https://www.health.qld.gov.au/mass/eapply>



# Questions and Feedback



Fill in the feedback form [here](#) to receive a personalised certificate of attendance



# Thank you

[MASS-Education@health.qld.gov.au](mailto:MASS-Education@health.qld.gov.au)

[MASS-Equipment@health.qld.gov.au](mailto:MASS-Equipment@health.qld.gov.au)

[MASS-Equipment-TSV@health.qld.gov.au](mailto:MASS-Equipment-TSV@health.qld.gov.au)

