

Tena Pads InstaDry

Continence Aids Product Information Sheet

Tena Pads InstaDry



Description of Product: Disposable Adhesive Shaped Pads are held in place by an adhesive strip that adheres to firm fitting underwear. Designed to instantly lock liquid away.

Available in:

Product Name	Product Code	Size	Colour	Total Capacity	Current MASS Distributor
Tena Pads InstaDry Standard Length	760852*	12.1x31.2cm	white	585mL	Brightsky
Tena InstaDry Long Length	760860	12.4x33.8cm	white	600mL	Brightsky

*Tena Sample website has Tena Pads Instadry Standard length code as 760851 (10 pack). This is the same MASS product as 760852 (20 pack).



© State of Queensland (Queensland Health) 2024
creativecommons.org/licenses/by/3.0/au

For more information contact:

Medical Aids Subsidy Scheme
 PO Box 281
 Cannon Hill QLD 4170 Australia
 tel (07) 3136 3665

MASS-ContinenceAids@health.qld.gov.au
health.qld.gov.au/mass