

# Hollister Inview Leg Bag

## Continance Aids Product Information Sheet

### Hollister Inview Leg Bag - Sterile



**Description of Product:** Hollister Inview Leg Drainage Bags attach directly to a sheath, indwelling catheter or night drainage bag. The bag features a non-return valve and fabric backing.

NB: Leg straps are supplied with these bags.

**Available in:**

Product Name	Product Code	Outlet Tap	Tubing Length	Total Capacity	Current MASS Distributor
Hollister Inview Leg Bag	9621	Lever Tap	10cm	500mL	Independence
Hollister Inview Leg Bag	9624	Lever Tap	50cm	500mL	Independence
Hollister Inview Leg Bag	9631	Lever Tap	10cm	800mL	Independence
Hollister Inview Leg Bag	9632	Lever Tap	50cm	800mL	Independence



© State of Queensland (Queensland Health) 2024  
[creativecommons.org/licenses/by/3.0/au](https://creativecommons.org/licenses/by/3.0/au)

**For more information contact:**

Medical Aids Subsidy Scheme

PO Box 281

Cannon Hill QLD 4170 Australia

tel (07) 3136 3665

[MASS-ContinanceAids@health.qld.gov.au](mailto:MASS-ContinanceAids@health.qld.gov.au)

[health.qld.gov.au/mass](https://health.qld.gov.au/mass)