

# Hollister Incare

## Continance Aids Product Information Sheet

### Hollister Incare



**Description of Product:** Hollister Incare Night Drainage Bags are sterile rectangular shaped bags that attach directly to a sheath, indwelling catheter or leg bag. This bag features a non-return valve.

**Available in:**

Product Name	Product Code	Outlet Tap	Tubing Length	Hanging Method	Total Capacity	Current MASS Distributor
Hollister Incare	9650	Lever Tap	120cm	No hanger included	2L	CH2



© State of Queensland (Queensland Health) 2024  
[creativecommons.org/licenses/by/3.0/au](https://creativecommons.org/licenses/by/3.0/au)

**For more information contact:**

Medical Aids Subsidy Scheme

PO Box 281

Cannon Hill QLD 4170 Australia

tel (07) 3136 3665

[MASS-ContinenceAids@health.qld.gov.au](mailto:MASS-ContinenceAids@health.qld.gov.au)

[health.qld.gov.au/mass](https://health.qld.gov.au/mass)